

GLOSEAL™

Fluorescent Dye with Sealant

- ▶ Pinpoint micro-leaks with our exclusive, **OEM-Grade** dye formula
- ▶ **2-in-1** premium dye with powerful sealant
- ▶ **Concentrated dye** ensures all leaks are found
- ▶ Sealant repairs hidden micro-leaks and stops deterioration from the inside
- ▶ **Perfect for preventative maintenance** – Scan the system to reveal new leaks
- ▶ **Works for you 24/7** – Continuous leak detection and ongoing sealing



PREMIUM DYE

+ POWERFUL SEALANT

1. Find & Fix
Immediate Leaks



2. Seal Difficult
to Find Leaks



3. Leak Detection &
Ongoing Sealing

THE POWER OF GLOSEAL™

OEM-GRADE
Fluorescent Dye

2-IN-1
Solution





COMPATIBLE
with All Refrigerants

POWERFUL & SAFE
Non-Polymer Sealant,
Not Activated by Air or Moisture

24/7
Performance

MULTIPLE
Delivery Methods

How GLO Seal™ Works



1. Find & Fix Immediate Leaks

Scan the system with a UV Light to find all leak areas glow brightly with GLO Seal™ concentrated dye. Make an external repair to ensure all large leaks are fixed. GLO Seal™ can remain safely in the system.



2. Seal Difficult to Find Leaks

GLO Seal™ pinpoints the location of all internal micro-leaks, forming a repair from the inside.



3. 24/7 Leak Detect + Ongoing Sealant

Scan the system periodically to find new leaks. Small, hard-to-find leaks hurt system performance and will eventually lead to total compressor burnout. Have peace of mind as GLO Seal™ works for you 24/7 for continual leak detection and ongoing sealing!

EZ-Ject™ Injection Method



GLO Seal™ Fluorescent Dye + Sealant EZ-Ject™ Dye Injection Kit

SPE-EZDSK-CS

- ▶ SPE-EZDS (2) 0.5 oz (15 ml) EZ-Ject™ cartridges prefilled with GLO Seal™ Fluorescent Dye + Sealant
- ▶ RP-EZ-50 EZ-Ject™ dye injector assembly
- ▶ RP-EZ-25 8 in (20 cm) EZ-Ject™ hose assembly
- ▶ RP-AD-410A R-410A adapter/purge fitting
- ▶ Service labels

Treats a **TOTAL** of up to:
▶ **7 tons** (24.6 kW) of cooling



Reusable Injector Assembly:
Save money on unnecessary hardware with each application!




GLO Seal™ Fluorescent Dye + Sealant EZ-Ject™ Cartridge

SPE-EZDS-CS

- ▶ 0.5 oz (15 ml) EZ-Ject™ cartridge prefilled with GLO Seal™ Fluorescent Dye + Sealant

Treats up to:
▶ **3.5 tons** (12.3 kW) of cooling

Connection:



Stick Capsule Injection Method




GLO Seal™ Fluorescent Dye + Sealant Stick Capsule

SPE-GSDS-CS

- ▶ 0.5 oz (15 ml) stick capsule prefilled with GLO Seal™ Fluorescent Dye + Sealant
- ▶ Requires bleed valves for injection

Treats up to:
▶ **3.5 tons** (12.3 kW) of cooling

Connection:



Bleed Valve

RP-BV-50

- ▶ Requires two (not included)



265 Spagnoli Rd, #100, Melville, NY 11747
(516) 333.4840 | Orders@Spectroline.com
CustomerService@Spectroline.com

www.Spectroline.com



G-41, Varalaxmi Complex,
MG Road, Secunderabad 500003
wrc.kaushal@gmail.com
9701525128
www.wrvcvac.com

SCAN QR
CODE
TO VIEW
VIDEO



HVAC/R
01/22 A19254-7
PRINTED IN USA